

KeraScan Screen-Positive Referral Report

Report identifier: 8C38E6E094065BB7 | Generated: 2026-08-21 15:26 UTC

Screening information

Field	Value
Anonymous child/study ID	SYNTHETIC-EXAMPLE-001
Screening date and time	2026-08-21T09:30:00Z
Operator ID	SYNTHETIC-OPERATOR
Device/configuration version	SYNTHETIC-DEVICE / synthetic-example-no-model
Protocol version	kerascan-school-screening-provisional-1
Model version	synthetic-example-no-model
Software version	phase4-school-screening
Image-analysis provenance hash	OD: 5cc3971118dc8a674a47812f98b2d8b52b0d58d53c7367b97bfed18a60b5477e

Final outcome

Child-level result	Referral priority	Affected eye
SCREEN-POSITIVE	PRIORITY_1	OD

Why the child screened positive

Eye	Screening domain	Observed value/result	Study threshold	Reason
OD	KeraScan image	Suspicious	—	Image pattern classified as suspicious
OD	K2	47.20 D	>46.8 D	Above study threshold


The approved image-screening model classified the KeraScan pattern as suspicious. Detailed feature attribution is not available for this model version.

Complete simplified measurements


Eye	K1	K2 / status	Pachymetry / status	Cylinder / magnitude status	Image status	Per-eye decision
OD	43.20 D	47.20 D / ABNORMAL	520 µm / WITHIN THRESHOLD	-0.50 D / WITHIN THRESHOLD	SUSPICIOUS	HIGH_RISK_SCREEN_POSITIVE
OS	42.80 D	44.10 D / WITHIN THRESHOLD	530 µm / WITHIN THRESHOLD	-0.25 D / WITHIN THRESHOLD	NORMAL_LIKE	SCREEN_NEGATIVE

KeraScan analysis images

OD analysis




OD: Cropped de-identified corneal/Placido ROI



OD: ROI with refined centre



OD: Tracked-ring overlay



OD: Directional-spacing visualization

Image-quality status	Tracking confidence	Geometry validation
ACCEPTABLE	0.91	PASS

KeraScan is an experimental initial screening aid. This report does not diagnose or exclude keratoconus. A screen-positive result indicates that further ophthalmic evaluation is recommended outside this application.